

**POWER**  
**SOCCER**  
**SHOP.COM**



A DIVISION OF **NEMI**

# *Quickie P200*

*Power Chair*

## *Steel Soccer Guard Mount*


*Installation Instructions*

PN 1628



POWER SOCCER SHOP INC.  
11841 243rd AVE NW  
ZIMMERMAN, MN 55398  
763-856-2044

# *Kit Contents*



Your new Quickie P200 soccer guard mount kit should contain the following items:

- (2) Main Guard Mounts
- (1) Guard Mount Crossmember
- (2) Ball Guards
- (2) Straps
- (1) Bolt Kit



Attach the right & left side guard mounts to the chair using (2) straps and (4) of the 2-1/4" bolts. Do not tighten bolts





*Crossmember Mounting Holes*

*Chair Crossmember*

*Guard Mount Crossmember*

Your new guard mount is designed to accommodate a wide variety of chair widths. You may or may not use all (8) bolts, depending on your chair configuration.

Position the guard mounts on the chair crossmember as such that the maximum number of bolts can be used to attach the guard mount crossmember to the guard mounts. Install the bolts, but do not tighten at this point.

Once the proper width is set and all bolts are in place, tighten all bolts *securely*.



Using the button head bolts, install the ball guards.